

Microscope Medical

Medical Microscope with Mechanical stage
With Optics with 100X Oil immersion
AB-05A

Huygenian Eye Piece	Achromatic Objective	Magnification
10X	10 X,45x	100 X,150x, 450X
15X	100 X oil immersion	675X,1000x,1500x

Viewing Head : Monocular tube
Body die cast : can be inclined body power coated
Nosepiece : Triple revolving Upright
Mechanical Stage : 120x120mm with stage mechanical
Focusing Knobs : Coarse and Fine separate
Condenser : Fixed Abbe NA 1.25 with Iris Diaphragm
Illumination : Plano-concave mirror
Carrying Case : Styrofoam Thermocole
Supplied with : Manual & Dust cove

Cat No	Type	Range	Optics
AB-05-A	Medical Microscope with mechanical stage	Iris diaphragm condenser lens focusable	Eye piece 10x,15x Objective 10x,45x without 100x Oil lens
AB-05-B	Medical Microscope with mechanical stage	Iris diaphragm condenser lens focusable	Eye piece 10x,15x Objective 10x,45x with 100x Oil lens



Abron Exports
www.abronexports.com

Optional Digital Camera Eye for Microscope

Cat No	with software for laptop	Mega Pixel
AB-14-A3	Digital Camera Eye for laptop with software	3 MP
	Resolution of display	2048*1536
AB-14-B5	Digital Camera Eye for laptop with software	5 MP
	Resolution of display	2560*1944
AB-14-C1	Digital Camera Eye for laptop with software	1.3 MP
	Resolution of display	1280*1024
AB-14-D10	Digital Camera Eye for laptop with software	10 MP
	Sensor & Pixel - 1/2.3" (6.44x4.616), 1.67x1.67um	
AB-14-E14	Digital Camera Eye for laptop with software	14 MP
	Sensor & Pixel 1/2.3"(1.4x1.4), upto 6.2-53fps,	

Microscope Lamp Abron

Cat No	Type	Range	Power
AB-57-B	Microscope LED Light sub stage	2x1.5V Cell	
AB-57-C	Microscope LED Light sub stage	Variable 2x1.5V Cell	
AB-57-D	Microscope LED Light sub stage	variable	220v
AB-59-B	Microscope Cleaning Kit	Brush, bellow, tissue, solu	
AB-60-A	Micrometer fitted in Eye Piece	1cm=	100 Part

Digi Microscope Medical

Medical Microscope Coaxial Mechanical stage
With Optics with 100X Oil immersion
AB-05B / Achromatic / Din Optics

Huygenian Eye Piece	Achromatic Objective	Magnification
10X	10 X,45x	100 X,150x, 450X
15X	100 X oil Immersion	675X,1000x,1500x

Viewing Head : Monocular tube
Body die cast : can be inclined body power coated
Nosepiece : Triple revolving Upright
Mechanical Stage : 120x120mm with coaxial stage
Focusing Knobs : Coarse and Fine separate
Condenser : Fixed Abbe NA 1.25 with Iris Diaphragm
Illumination : Plano-concave mirror
Carrying Case : Styrofoam Thermocole
Supplied with : Manual & Dust cover

Cat No	Type	Range	Optics
AB-05-C	Medical Microscope	Coaxial mechanical stage	Iris diaphragm condenser lens focusable Eye piece 10x,15x Objective 10x,45x 100x Oil lens
AB-06-E	Medical Microscope	LED Light Coaxial mechanical stage	Iris diaphragm condenser lens focusable Eye piece 10x,15x Objective 10x,45x 100x Oil lens
AB-06-G	Medical Microscope	LED Light 360 Tube + stage	Eye piece 10x,15x Objective 10x,45x 100x Oil lens

With special Din Optics

AB-06-D Medical Microscope Coaxial mechanical stage
Iris diaphragm condenser lens focusable

Din Optics Eye piece 10x,15x Objective 10x,45x 100x Oil lens



Abron Exports
www.abronexports.com

Microscope with Digital Camera Eye

Cat No	with software for laptop	Mega Pixel
AB-06-D1	Microscope Medical + Digital camera software	1.3MP
	Eye piece 10x,15x Objective 10x,45x,100x Oil lens	
AB-06-D3	Microscope Medical + Digital camera software	3 MP
	Eye piece 10x,15x Objective 10x,45x,100x Oil lens	
AB-06-D5	Microscope Medical + Digital camera software	5 MP
	Eye piece 10x,15x Objective 10x,45x,100x Oil lens	

Microscope Spares Abron

Cat No	Type	Long	Parts
AB-60-A	Micrometer fitted in Eye Piece	1cm=	100 Part
AB-66-A	Stage Micrometer Slide with	1cm=	100 Part
AB-66-B	Stage Micrometer Slide with	1mm=	100 Part
AB-65-A	Pointer Eye Piece	10X	
AB-64-A	Ocular Micrometer Disc with	1mm=	100 Part